





Your Health Matters

Use of Insulin Pen / 胰島素筆的使用方法

1. Remove protective tab / 取下筆尖保護套




2. Screw needle onto pen / 把針擰緊在筆上




- Remove outer and inner plastic needle caps
取下塑料針的內，外帽

3. Do air shot / 做空氣注射





- Dial to 2 / 撥到2
- Press button / 按下按鈕
- Repeat until you see insulin leak out of needle
重復，直到看見胰島素洩漏出針頭。
- Shake off excess insulin from tip
抖落尖頭多餘的胰島素

4. Dial the insulin dose / 撥到胰島素的劑量




5. Give the insulin / 給予胰島素



- Insert needle into skin at 90 degree angle
以90度直角把針插入皮膚
- Press button to deliver insulin
按下按鈕，注入胰島素
- Make sure dial returns to 0
確定轉盤回到0
- Count to 10 before removing needle
數到10，取針

6. Remove needle from pen / 把針從筆上取下



- Never store a needle on the pen as it will draw air into the pen
切勿把針留在筆上，以免把空氣抽入筆